

## 3 Foot Vertical lifting frame



The Carers bed mechanism can be inserted into any of our adjustable bed surrounds. It is especially good if the person has carers and needs to be turned in bed at a comfortable height but also to help the person stand up out of bed whilst sitting on the edge of the bed as it raises them to a standing position like a riser chair, you must sit in the middle section when raising the bed. The mechanism can raise up to 18 inches flat for turning the patient or raised whilst the bed is profiling (head and / or foot elevated).

The minimum starting height can be 21.5 inches from the floor to top of mattress, the standard starting height is 23 inches. (please see height adjustments below, under accessories).

The single vertical lifting frame can be put into any of our adjustable bed surrounds, so choose the bed you like (from the home page) and add the cost of this single vertical lifting frame to the cost of any bed (from dropdown mechanisms tab (7th tab down) on any bed page) to get the combination you want.

### Description

- 5 part adjustable profiling with integral raise and lower facility
- 4 motors (2 for profiling and 2 for lift)
- Vertical rise up to 18 inches from starting point
- Twin On-Line fittings
- Weight adjuster zone
- Centre belt
- On-Line shoes
- Load bearing points bolted and sleeved
- **Special order**

Options (See Bottom of Page)

- Massage units

**Dual beds** have independant controls and independant mattresses. You will generally have a vertical lifting mechanism on one side of a dual bed, so if you are purchasing a 5ft dual bed you will only need one 2ft 6 vertical lifting mechanism.

**Double beds** have one mattresses, two vertical lifting mechanisms and a special linking cabel.

**THE PRICE OF THE VERTICAL LIFTING FRAME WILL BE ADDED TO THE COST OF YOUR BED.**

This price includes delivery & set up by trained engineers.

€995.00

[Vendor Information](#)